





# Product Comparison

the medcom group, ltd

|  | ThermX                                                                              | Game Ready                                                                          | Vascutherm 5                                                                         | Nice 1                                                                                | Breg Polar Care Wave      | Breg Polar Cube                       | Cryo Cuff (Gravity) | Cryo Cuff (IC)     | Dr. Active CCT               | DonJoy Iceman           | ARTEK Sport   |
|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------|---------------------------------------|---------------------|--------------------|------------------------------|-------------------------|---------------|
| Rental Option                                                                     | ✓                                                                                   | ✓                                                                                   | ✓                                                                                    | ✓                                                                                     | ✗                         | ✗                                     | ✗                   | ✗                  | ✗                            | ✗                       | ✗             |
| RX Needed                                                                         | ✓                                                                                   | ✓                                                                                   | ✓                                                                                    | ✓                                                                                     | ✓                         | ✓                                     | ✗                   | ✓                  | ✗                            | ✓                       | ✓             |
| Iceless                                                                           | ✓                                                                                   | ✗                                                                                   | ✓                                                                                    | ✓                                                                                     | ✗                         | ✗                                     | ✗                   | ✗                  | ✓                            | ✗                       | ✗             |
| Cold Temperature Range                                                            | 38°F - 55°F                                                                         | 35°F - 50°F                                                                         | 43° - 49°F                                                                           | 43°F - 59°F                                                                           | 45°F-50°F                 | 35°F - 50°F                           | 35°F - 50°F         | 35°F - 50°F        | N/A                          | 35°F - 50°F             | 35°F - 50°F   |
| Wraps Available                                                                   | Knee, Hip, Shoulder, Lumbar, Universal, Ankle                                       | Knee, Articulated Knee, Shoulder, Hip, Ankle                                        | Knee, Shoulder, Hip, Elbow, Foot/Ankle, Hand, Lumbar, Cervical                       | Knee, Articulated Knee, Shoulder, Hip, Ankle, Lumbar, Butterfly                       | Knee, Shoulder, Universal | Knee, Shoulder, Hip, Ankle, Butterfly | Knee, Ankle, Elbow  | Knee, Ankle, Elbow | Knee, Shoulder, Ankle, Wrist | Shoulder, Universal     | Universal     |
| Weight (empty)                                                                    | 14 lbs                                                                              | 8 lbs                                                                               | 18 lbs                                                                               | 10 lbs                                                                                | 5 lbs                     | 8 lbs                                 | 5 lbs               | 5 lbs              | 2 lbs                        | 5 lbs                   | 4 lbs         |
| Dimensions                                                                        | 10" x 9" x 15"                                                                      | 5.5" x 16" x 8.5"                                                                   | 6.6" x 12" x 13.2"                                                                   | 8" x 7.64" x 8"                                                                       | 12" x 12"                 | 12" x 12"                             | N/A                 | N/A                | N/A                          | 12.75" x 12.5" x 12.25" | 10" x 8" x 8" |
| Compression Option                                                                | ✓                                                                                   | ✓                                                                                   | ✓                                                                                    | ✓                                                                                     | ✓                         | ✗                                     | ✓                   | ✓                  | ✓                            | ✗                       | ✗             |
| Compression Range                                                                 | 0-75 mmHG                                                                           | 0-75 mmHG                                                                           | 0-50 mmHG                                                                            | 0-75 mmHg                                                                             | 0-50 mmHg                 | N/A                                   | N/A                 | N/A                | N/A                          | N/A                     | N/A           |
| Sound Level Visual (Average highest level)                                        |  |  |  |  | N/A                       | N/A                                   | N/A                 | N/A                | N/A                          | N/A                     | N/A           |
| Sound Level dBa (Average highest level)                                           | 56 dBa                                                                              | 65 dBa                                                                              | 57 dBa                                                                               | 75 dBa                                                                                | N/A                       | N/A                                   | N/A                 | N/A                | N/A                          | N/A                     | N/A           |
| Heat Therapy                                                                      | ✓                                                                                   | ✗                                                                                   | ✓                                                                                    | ✗                                                                                     | ✗                         | ✗                                     | ✗                   | ✗                  | ✗                            | ✗                       | ✗             |
| Heat Temperature                                                                  | 105°                                                                                | N/A                                                                                 | 105°                                                                                 | N/A                                                                                   | N/A                       | N/A                                   | N/A                 | N/A                | N/A                          | N/A                     | N/A           |
| Contrast Therapy                                                                  | ✗                                                                                   | ✗                                                                                   | ✓                                                                                    | ✗                                                                                     | ✗                         | ✗                                     | ✗                   | ✗                  | ✗                            | ✗                       | ✗             |
| Built in Timer                                                                    | ✓                                                                                   | ✓                                                                                   | ✓                                                                                    | ✓                                                                                     | ✗                         | ✗                                     | ✗                   | ✗                  | ✗                            | ✗                       | ✓             |
| Touch Screen                                                                      | ✓                                                                                   | ✗                                                                                   | ✓                                                                                    | ✓                                                                                     | ✗                         | ✗                                     | ✗                   | ✗                  | ✗                            | ✗                       | ✓             |
| DVT Sleeve Compatible                                                             | ✓                                                                                   | ✗                                                                                   | ✓                                                                                    | ✗                                                                                     | ✗                         | ✗                                     | ✗                   | ✗                  | ✗                            | ✗                       | ✗             |
| Bilateral Wrap Compatible                                                         | ✗                                                                                   | ✓                                                                                   | ✗                                                                                    | ✗                                                                                     | ✗                         | ✓                                     | ✗                   | ✗                  | ✗                            | ✗                       | ✗             |